


Agrosta® Winterwood

Durofel Firmness and size measurement of soft fruits

The Agrosta®Wonderfast / Durofel version was designed in 2020 to provide a highly accurate and fully automated solution for testing the firmness of Cherries, Blueberries, Tomatoes, Cranberries, and Grapes.



Thanks !

Many thanks for having acquired an Agrosta instrument
Your package contains :

- The Agrosta®Wonderfast / Durofel machine
- 2 turntables (Depending on your request, either two identical turntables or two different ones)
- A user manual with a certificate of conformity
- A 110V / 220V power supply
- A USB stick containing the software
- A USB cable

We highly recommend reading the user manual, as it provides clear and concise instructions to help you maximize the use of the device.

Product Overview

Design & Manufacturing : The **Agrosta®Winterfast / Durofel** version is entirely designed and assembled in France.

- The **motherboards** come from the **USA**.
- The **shield boards** are manufactured in **Hong Kong and the USA**.

Handling & Care: The **Agrosta®Winterfast / Durofel** version is **NOT waterproof**. It is a **precision instrument**, so please handle it with care and avoid **dropping or knocking** it.

Warranty: The **Agrosta®Winterfast / Durofel** version comes with a **two-year guarantee**.

Measurement Specifications :

- **Minimum Durofel:** 7%
- **Maximum Durofel:** 100% (*Shore A Scale*)
- **Resolution:** ±1%
- **Size Measurement:** Provided in **mm** (*The user must enter the initial size between the sensor head and the bottom of the cup of the turntable*)

Compatibility : The software is compatible with Windows XP, Vista, 7, 8, 10 and 11

Fruit Measurement Impact : Fruits like **Blueberries, Blackberries, Strawberries, and Raspberries** are destroyed during measurement. This **does not affect the results**, as the machine records **only the maximum pressure**.

AGROSTA®Winterwood – Step-by-Step Guide

1. Install the Driver

- Do not connect your machine yet.
- Insert the **USB stick** into your computer.
- In **2023**, there are **two drivers** to install:
 1. **Install DRIVER1 first.**
 2. **Then install DRIVER2.**

 Agrosta Winterwood.EXE	29/04/2021 09:...	Application	33 769...
 Agrosta_Driver.EXE	24/01/2017 01:...	Application	238 Ko

2. Connect the USB Cable

- Plug the **USB cable** between the **instrument** and your **computer**.

3. Device Recognition

- **Wait a few seconds** until the device is recognized.
- The **driver will automatically link to the device**.

4. Install the Software

- Open the **USB stick** and **double-click** on **"INSTALL"**.
- **Follow the setup procedure** to complete the installation.


 Agrosta Winterwood.EXE	29/04/2021 09:...	Application	33 769...
 Agrosta_Driver.EXE	24/01/2017 01:...	Application	238 Ko

5. Connect the Power Plug

- Plug the **power supply** into the **machine**.

6. Start the Software

- Once installed, the **software will start automatically**.
- The **connection with the machine is automatic**.



7. Software Usage & Testing Procedure

The software comes with a **light version of Excel**, allowing you to:

- Open any **Excel file** or use the **blank sheet** that appears at startup.
- Click on the cell where you want the **data to start displaying**.
- Click on **"START NEW SERIE"**: This will add the **date, time, and column headers**.
- You can **click on any other cell** to add extra information (*e.g., sugar content, variety...*) or start another batch in a different part of the sheet.
- At the end of your tests, you can **save the Excel file**.


8. Barcode Reader Support

- You can use a **barcode reader** to scan references.
 - **Select the target cell**, then scan the barcode to insert its reference.
-

9. Launching a Cycle

1. Define the number of tests

- Enter the **number of fruits** placed on the table (between **1 and 32**) in the corresponding field.



The image shows a screenshot of a software interface with two input fields. The first field is labeled "INITIAL SIZE mm" and contains the value "34". The second field is labeled "NUMBER OF TESTS" and contains the value "2".


-
- Example for **blueberries**:
 - If the table is **full of fruits**, enter **32**.
 - For **continuous testing** (*removing and replacing berries during the cycle*), enter a **higher value** like **1000**.
 -

2. If you need to stop the cycle

- You can **stop anytime** by clicking on **"STOP CYCLE"**.
- If the machine does not stop properly, unplug the **USB and power cable**, then reconnect and restart the software.

3. Position the turntable

- Switch between **two cups** as shown in the reference photo.
- The machine **always starts** with the **cup on the right side**.



4. Start the cycle

- Click on **"LAUNCH CYCLE"**.
- The machine will begin testing based on the selected parameters.

WINTERWOOD by AGROSTA

START NEW SERIE
Click before measurement

STATISTICS
Click after measurement

FIRMNESS

SIZE

IN MM ADJUSTED WITH INITIAL SIZE

Home Insert Formulas Layout Display

New Open Save Print Paste Copy

Verdana 10 A⁺ A⁻

Automatic line wrap

B I U Merge and center

Format Cells Editing

Document Clipboard Font Alignment Number Styles

	A	B	C	D	E	F	G	H	I	J	K	L	M	N
1														
2														
3														
4														
5														
6														
7														
8														
9														
10														
11														
12														
13														
14														
15														
16														
17														
18														
19														
20														
21														
22														
23														
24														
25														
26														
27														

Worksheet1

INITIAL SIZE mm

NUMBER OF TESTS

LAUNCH CYCLE

10. Adjusting the Size Measurement


To ensure accurate fruit size measurement, you need to **enter the initial size** in millimeters.

Setting the Initial Size


- This value corresponds to the **distance between the sensor head and the middle of the cup height** of the turntable you are using.
- Enter this value in the "**Initial size**" field in the software.

Calibration & Adjustment

- Perform **test measurements** with fruits of a **known size**.
- If needed, adjust the initial value:
 - **Adding 1 mm to the initial size will add 1 mm to the measured fruit size.**
 - Adjust accordingly until the measurements match the actual fruit sizes.
 -



INITIAL SIZE mm	34
NUMBER OF TESTS	2



Differences in Measurement Values Between the Handheld Agrosta®100 and the Agrosta®Winterfast / Durofel Version



Measurement Differences on Blueberries

The **difference in measurements** occurs **only on blueberries**, due to the **handling method**:

- **With handheld devices (Agrosta®100):**
 - The blueberries are held between **two fingers**, which applies **overpressure**, leading to **higher values**.
- **With the Agrosta®Winterfast / Durofel version:**
 - The blueberries are placed **freely inside the turntable cups**, resulting in **lower values**.
 - The difference is approximately **16%** based on experience.
 - To **compensate**, multiply the values obtained by **1.16**.
 - However, the variation depends on how users **hold the fruit** with handheld devices.

Effect of Fruit Breaking During Measurement

- **Mature blueberries may be crushed** by the machine during the test.
- This **does not affect the measurement**, because:
 - The machine **detects** when it touches the fruit.
 - It **records the maximum pressure every 100 milliseconds**.
 - If the fruit **breaks**, the pressure drops, but the **maximum pressure is already stored**.

Cleaning & Maintenance

- Remove and **clean the tip daily** to prevent fruit juice from **blocking the system**.
- **Be careful** when **screwing and unscrewing** the **plastic abutment** to avoid damage.



Minimum Fruit Height Requirement



- The fruit **must exceed the turntable level by at least 3 mm**.
- If the fruit is **too small**, the tray may **touch the endstop before completing the measurement**.

How the Sensor Moves

1. The **sensor detects the fruit**.
2. It must **move down approximately 3 mm** after touching the fruit.
3. If the sensor **touches the endstop** after only **1 or 2 mm**, the tray **goes back before finishing the measurement**.

"TOO SMALL" Error Message

- If the sensor **cannot move 3 mm**, the software displays **"TOO SMALL"**, indicating an **incomplete measurement**.
- Ensure that the fruit meets the **minimum height requirement** to avoid this error.

4 / 29 / 2020		15 : 36
Durofel %	Size	
37	20	
48	23	
TOO SMALL	11	
41	20	
40	16	
31	16	
43	15	
41	22	
41	21	
44	18	

Solution for Incomplete Measurements

If you encounter the **"TOO SMALL"** error or obtain **lower values**, the recommended solution is:

- **Use a turntable with smaller cups** to ensure proper support for small berries.
- When placing **small berries**, **block each fruit** in the cup to prevent it from **moving down during measurement**, as this leads to **inaccurate and lower values**.



Ensuring Proper Blueberry Positioning

- When placing **blueberries in the cups**, make sure to **LOCK each berry carefully**.
 - If a berry moves **even slightly downward** during measurement, the **Durofel Index will be lower than expected**.
 - Example:
 - A **0.25 mm movement** results in a **10% decrease** in the Durofel Index.
 - This is why the device is supplied with **turntables of different cup sizes**, allowing proper adaptation for various fruit sizes.
-

Cleaning & Maintenance

Preventing Fruit Juice Contamination

To maintain the precision of the machine, always keep it **free of fruit juice**.

Cleaning Procedure

1. **Unscrew the abutment** and remove the tip.
2. **Clean the tip and tip chamber** thoroughly with **hot water**.

Using the Tip Protection (Protective Covering)

- The **protective covering** must be placed **before starting the machine** (see *reference photo*).
- It should **neither be too loose nor too tight**.
- **For Blueberries**, using the protection is **mandatory**.
- **For other fruits**, it is **highly recommended**.



Sensor Protection & Juice Management



- The **sensor is highly sensitive to fruit juice.**
- You have two cleaning options:
 1. **Use a finger cot (recommended)** – A **standard disposable finger cot** found in most countries.
 2. **Manually clean after each juice spill** – This involves:
 - **Unscrewing the abutment.**
 - **Cleaning the tip, chamber, and abutment each time a fruit is broken and juice flows.**
- **Important:** If using a finger cot, replace it **immediately if it is damaged** to ensure proper protection.



Calibration Modification (For Experts Only)

⚠ **Warning:** Any modification in the calibration settings is **permanent**, and the machine will lose its previous calibration.



Calibration Check Procedure


If you follow the procedure below and obtain a **Durofel value between 96% and 100%**, your machine is correctly calibrated.

Why is this test reliable?

The **Winterwood machine** follows the **Durofel Standard**, which itself is based on the **Shore A hardness scale**.

- In the **Shore A Standard**, the tip (in green and blue) is linked to a **calibrated spring** that applies a **precise load**.
- The **tip retracts** depending on the **hardness of the sample**.

Durometer hardness test



How Calibration Works in the Winterwood Sensor Head

- The **calibrated spring** is attached to a **calibrated load cell**.
- The **load cell measures the applied pressure**, which is interpreted as a **retraction length**.
- **100% Durofel** corresponds to **806 grams of pressure** when the tip is **fully retracted**.

Potential Calibration Failures

1. Calibrated Spring Failure:

- If the spring becomes **too soft**, the tip **will not reach 806 grams of pressure** at full retraction.
- The machine will **not exceed 96% Durofel**.

2. Calibrated Load Cell Failure:

- If the **load cell provides incorrect values**, the pressure will never reach the correct level.


✓ The only way to get a value between 96% and 100% with a faulty calibration would be if **both components failed in opposite ways and compensated each other**.

✗ This scenario is **highly unlikely** because **wear typically causes both components to drift in the same direction**.

How to Check Calibration?

Step-by-Step Calibration Test

1. **Prepare a flat, hard surface** (or use the calibration template if provided).
2. **Remove the tip protection** (if installed).
3. **Remove the turntable.**
4. **Open the software** and **connect the machine**.
5. **Start a measurement cycle.**
6. **During the sensor tray movement downward:**
 - Place the **flat piece** in front of the sensor.
 - **Press firmly** and hold for **one second**.
7. **Remove the flat piece** before the sensor tray reaches the bottom of the machine.
8. **Read the measured value:**
 - If it is **between 96% and 100% Durofel**, the calibration is **correct**.
 - If it is **below 96%**, the machine needs **calibration adjustment** or **maintenance**.



Optional Colorimeter (Winterwood Model)

The color measurement is performed **from under the table**.

The screenshot displays the WONDERFAST software interface. The main window shows a data table with columns for color measurements (Purple, Blue, Green, Yellow, Orange, Red, Ctrial Cherry) and rows for individual measurements and averages. The right-hand side features a control panel with various settings and a 'LAUNCH CYCLE' button.

	Durofel %	Size	Purple	Blue	Green	Yellow	Orange	Red	Ctrial Cherry
7	76	21	25	10	40	27	66	68	1
8	68	25	24	10	40	26	66	67	1
9	69	28	20	6	39	21	60	57	2
10	64	23	21	6	38	22	59	57	2
11	67	34	23	5	38	19	49	44	3
12	74	24	24	7	39	20	50	45	3
13	66	26	20	4	38	19	47	30	4
14	58	30	18	4	38	19	47	30	4
15	75	28	19	4	37	18	46	11	5
16	70	29	18	4	38	19	47	11	5
17	64	34	22	5	38	18	45	6	6
18	71	20	18	4	35	18	41	6	6
19	70	19	16	3	33	16	38	4	7
20	68	19	18	3	37	18	42	3	7
21	AVERAGE / STDEV								
22	68.57								
23	4.8								

Color Measurement Details

- The device measures **six different colors**.
- These values are **relative** and follow a scale where:
 - **100 = Pure White** (Maximum value)
 - **0 = Pure Black** (Minimum value)
- The measured **wavelengths** correspond to:
 - **450 nm**
 - **500 nm**
 - **550 nm**
 - **570 nm**
 - **600 nm**
 - **650 nm**



Excel Functions in the Software

The software includes a **built-in light version of Excel**, allowing you to use all **standard Excel functions** such as:

- **=AVERAGE()** → Calculate the average of a dataset.
- **=STDEV()** → Calculate the standard deviation.
- **Many other functions** are available in the **"Expression"** tab.

Quickly Access Batch Statistics

- When you have **finished a batch**, you can obtain **quick statistics** by clicking on the **"STATISTICS"** button, located at the **top center of the window**.
-  **This feature was introduced in 2021.**



Congratulations! 🎉 You are now ready to fully utilize the **AGROSTA® Winterwood**.

If needed, refer to the user manual for additional guidance and explore all the features offered by this device.

Happy measuring! 🍇🍒