

# Agrosta® Wonderfast

Digital Firmness and size measurement of fruits

---

**The Agrosta® Wonderfast was first designed in 2020 and improved in 2024 to provide a highly accurate and fully automated solution for testing any kind of fruit using penetrometry *and* Deflection.**



## Thanks !

Many thanks for having acquired an Agrosta instrument  
Your package contains :

- The Agrosta®Wonderfast
- 2 turntables (Depending on your request, either two identical turntables or two different ones)
- 2 Tips : For penetrometry standard or according to client request
- A user manual with a certificate of conformity
- A 110V / 220V power supply
- A USB stick containing the software
- A USB cable

**We highly recommend reading the user manual, as it provides clear and concise instructions to help you maximize the use of the device.**

## Product Overview

**Design & Manufacturing :** The **Agrosta®Wonderfast** is entirely designed and assembled in France.

- The **motherboards** come from the **USA**.
- The **shield boards** are manufactured in **Hong Kong and the USA**.

**Handling & Care:** The **Agrosta®Wonderfast** is **NOT waterproof**. It is a **precision instrument**, so please handle it with care and avoid **dropping or knocking** it.

**Warranty:** The **Agrosta®Wonderfast** comes with a **two-year guarantee**.

### Measurement Specifications :

- **Minimum Pressure : 100 grams**
- **Maximum Pressure: 10 000 grams**
- **Resolution: ±1g**
- **Size Measurement:** Provided in **mm** (*The user must enter the initial size between the sensor head and the bottom of the cup of the turntable*)

**Compatibility :** The software is compatible with Windows XP, Vista, 7, 8, 10 and 11

# AGROSTA®Wonderfast– Step-by-Step Guide

## 1. Install the Driver

- Do not connect your machine yet.
- Insert the **USB stick** into your computer.
- In **2023**, there are **two drivers** to install:
  1. **Install DRIVER1 first.**
  2. **Then install DRIVER2.**

 Agrosta Winterwood.EXE	29/04/2021 09:...	Application	33 769...
 Agrosta_Driver.EXE	24/01/2017 01:...	Application	238 Ko

## 2. Connect the USB Cable

- Plug the **USB cable** between the **instrument** and your **computer**.

## 3. Device Recognition

- **Wait a few seconds** until the device is recognized.
- The **driver will automatically link to the device**.

## 4. Install the Software

- Open the **USB stick** and **double-click** on **"INSTALL"**.
- **Follow the setup procedure** to complete the installation.


 Agrosta Winterwood.EXE	29/04/2021 09:...	Application	33 769...
 Agrosta_Driver.EXE	24/01/2017 01:...	Application	238 Ko

## 5. Connect the Power Plug

- Plug the **power supply** into the **machine**.

## 6. Start the Software

- Once installed, the **software will start automatically**.
- The **connection with the machine is automatic**.



## 7. Software Usage & Testing Procedure

The software comes with a **light version of Excel**, allowing you to:

- Open any **Excel file** or use the **blank sheet** that appears at startup.
- Click on the cell where you want the **data to start displaying**.
- Click on **"START NEW SERIE"**: This will add the **date, time, and column headers**.
- You can **click on any other cell** to add extra information (*e.g., sugar content, variety...*) or start another batch in a different part of the sheet.
- At the end of your tests, you can **save the Excel file**.


## 8. Barcode Reader Support

- You can use a **barcode reader** to scan references.
  - **Select the target cell**, then scan the barcode to insert its reference.
- 

## 9. Launching a Cycle

### 1. Define the number of tests

- Enter the **number of fruits** placed on the table (between **1 and 32**) in the corresponding field.



The image shows a screenshot of a software interface with two input fields. The first field is labeled "INITIAL SIZE mm" and contains the number "34". The second field is labeled "NUMBER OF TESTS" and contains the number "2".


- 
- Example for **blueberries**:
  - If the table is **full of fruits**, enter **32**.
  - For **continuous testing** (*removing and replacing berries during the cycle*), enter a **higher value** like **1000**.
  -

### 2. If you need to stop the cycle

- You can **stop anytime** by clicking on **"STOP CYCLE"**.
- If the machine does not stop properly, unplug the **USB and power cable**, then reconnect and restart the software.


### 3. Position the turntable

- Switch between **two cups** as shown in the reference photo.
- The machine **always starts** with the **cup on the right side**.



#### **Warning**

- The **largest fruit** must be placed **first** at the start of the cycle.
  - This will **define the maximum size** for the machine's operation.
- In the following example, the **first fruit to be tested** is the **pomelo** (*the largest fruit in the batch*).



#### 4. Measurement Strategies

The **Agrosta® Wonderfast** allows you to choose between **two different measurement strategies**, depending on the type of fruit and the desired analysis method:

##### 1 Penetrometry Measurement (Fixed Travel / Variable Pressure)

- In this method, you **set the penetration depth (RUN/TRAVEL)** inside the fruit.
- The machine will then **measure the pressure required** to reach this depth.
- This method is commonly used for **penetrometry tests** to evaluate fruit firmness.
- **Fill in the fields highlighted in YELLOW** in the software.

##### 2 Deformation Measurement (Fixed Pressure / Variable Deformation)

- In this method, you **set a fixed pressure**, and the machine will measure the **deformation** of the fruit. This strategy is typically used for **soft fruits**. It allows measurements according to the **FIRMTECH standard**:
  - **Low pressure** corresponds to the **TRIGGER**.
  - **High pressure** corresponds to **PRESSURE GRAMS**.
- **Fill in the fields highlighted in ORANGE** in the software.

The screenshot displays the software interface for the Agrosta®100 BLUE. It features a menu bar (File, Edit, Expression, Layout, Display), a toolbar with various icons, and a main workspace divided into a data table and a control panel.

**Data Table:**

Measurement	Size	Purple %	Blue %	Green %	Yellow %	Orange %	Red %
1 / 3 / 2023	9 : 4						

**Control Panel:**

- START NEW SERIES** (Click before measurement)
- STATISTICS** (Click after measurement)
- RESULT** (IN GRAMS OR IN 0.1MM)
- SIZE** (IN MM ADJUSTED WITH INITIAL SIZE)
- INITIAL SIZE mm** (140)
- TRIGGER (GRAMS)** (100)
- NUMBER OF TESTS** (3)
- Fix PRESSURE & Measure TRAVEL** (Set cherry parameters, Set blueberry param.)
- Fix TRAVEL & Measure PRESSURE** (Set tomato parameters, Set apple parameters, Set cherry parameters, Set blueberry param., Set strawberry param.)
- SPEED IN 0.1 MM/S** (150)
- PRESSURE GRAMS** (500)
- SPEED IN 0.1 MM/S** (50)
- RUN IN 0.1 MM** (60)
- LAUNCH CYCLE (0.1mm)**
- LAUNCH CYCLE (Grams)**

Red arrows in the image point to the highlighted fields: yellow for 'INITIAL SIZE mm', 'TRIGGER (GRAMS)', and 'NUMBER OF TESTS'; orange for 'SPEED IN 0.1 MM/S', 'PRESSURE GRAMS', and 'RUN IN 0.1 MM'.

## 10. Adjusting the Size Measurement

To ensure accurate fruit size measurement, you need to **enter the initial size** in millimeters.

### Setting the Initial Size

- This value corresponds to the **distance between the sensor head and the middle of the cup height** of the turntable you are using.
- Enter this value in the **"Initial size"** field in the software.


### Calibration & Adjustment

- Perform **test measurements** with fruits of a **known size**.
- If needed, adjust the initial value:
  - **Adding 1 mm to the initial size will add 1 mm to the measured fruit size.**
  - Adjust accordingly until the measurements match the actual fruit sizes.
  -



INITIAL SIZE mm

NUMBER OF TESTS



## 11. Adjusting the Trigger, Speed, and Run (Stroke) Parameters (In RUN Adjustment Mode)

TRIGGER (GRAMS)	32
NUMBER OF TESTS	3
SPEED IN 0.1 MM/S	30
RUN IN 0.1 MM	45

The **RUN adjustment mode** allows you to fine-tune the key parameters that determine how the **Agrosta® Wonderfast** interacts with the fruit during testing.

### ◆ TRIGGER – Contact Pressure Threshold

- The **TRIGGER** defines the **pressure** at which the machine **detects that it has touched the fruit**.
- Once the **trigger pressure** is reached, the machine **begins the RUN** (penetration) as defined in **increments of 0.1 mm**.
- **Increasing the TRIGGER value** results in **deeper penetration** into the fruit.

**For fruits with an irregular surface** (e.g., raspberries, blackberries):

- Increase the **TRIGGER** value to **reduce the influence of surface irregularities**, ensuring more consistent measurements.

### ◆ SPEED & RUN (Stroke) – Post-Contact Parameters

- These two parameters control the **speed and depth** of penetration after the machine **touches the fruit**.
- Their **impact on measurement accuracy is very high**.

#### **SPEED Adjustments:**

- The maximum speed value is **200**, corresponding to **20 mm/second**.
- **If a value higher than 200 is entered**, the motor will block and stop moving (*this will not damage the machine*).
- There is **no minimum speed limit**.

#### **RUN (Stroke) Adjustments:**

- The **minimum RUN value is 20**, which corresponds to a **2 mm penetration depth**.
- If the **applied pressure reaches 13,000 grams**, the tray will **automatically return** without completing the full stroke.
- If you consistently get **exactly 13,000 grams**, you should:
  - **Reduce the RUN value** (*shorten the penetration stroke*).
  - **Use a smaller tip** for testing.

## Spring System

The **spring system** is designed to **hold the fruit in place** on the table during measurement.

- It helps **stabilize** the fruit to ensure **consistent** and **accurate** results.
- However, **it is not essential for the machine's operation.**

✔ **If you prefer, you can remove the spring system** without affecting the functionality of the device.



## Optional Colorimeter

The color measurement is performed from under the table.

	A	B	C	D	E	F	G	H	I	J	K
1	12 / 23 / 2022	15 : 51									
2	Measurement	Size	Purple %	Blue %	Green %	Yellow %	Orange %	Red %			
3	1325	132	0	0	-3	-1	4	0		No sample	
4	2509	107	95	91	93	99	102	101		White sheet	
5	2050	114	6	0	8	6	33	38		Cherry 1	
6	2125	94	5	0	10	4	36	38		Cherry 1	
7	2577	108	4	1	6	3	34	43		Cherry 2	
8	1758	95	3	0	6	1	23	32		Cherry 3	
9	14	110	7	1	10	6	41	50		Cherry 4	
10	12	107	8	3	9	7	40	50		Cherry 4	
11	7	54	53	52	41	57	46	57		Grey sample	
12	7	54	54	52	42	57	47	57		Grey sample	
13											
14											
15											
16											

The Agrosta® Wonderfast is capable of measuring six different colors, providing relative values based on the detected light reflection.

### Measurement Scale

- The maximum value is reached with "Pure White".
- The minimum value is assigned to "No Sample" (absence of fruit).

### Measured Wavelengths (in nanometers)

- 450 nm
- 500 nm
- 550 nm
- 570 nm
- 600 nm
- 650 nm

### Calibration Reference

- A pure white surface placed on the turntable will return a value of 1000 for each wavelength.

- Dominant colors in a fruit will correspond to the highest percentage values.
  - Example: Cherries produce high readings in the Red and Orange spectrum (*Wavelengths 570-650 nm*).

This system allows for accurate color classification and fruit quality analysis.



## Excel Functions in the Software

The software includes a **built-in light version of Excel**, allowing you to use all **standard Excel functions** such as:

- **=AVERAGE()** → Calculate the average of a dataset.
- **=STDEV()** → Calculate the standard deviation.
- **Many other functions** are available in the "Expression" tab.

## Quickly Access Batch Statistics

- When you have **finished a batch**, you can obtain **quick statistics** by clicking on the **"STATISTICS"** button, located at the **top center of the window**.
-  **This feature was introduced in 2021.**



**Congratulations!** 🎉 You are now ready to fully utilize the **AGROSTA® Wonderfast**.

If needed, refer to the user manual for additional guidance and explore all the features offered by this device.

**Happy measuring!** 🍇 🍒